




# The Fleischman Lab

Author	Sarah J. Morse		July 13, 2017																		
Title	CALR <sup>DEL</sup> and CALR <sup>INS</sup> Screening																				
Introduction	<p>Calreticulin (CALR) mutations are common in patients with essential thrombocythemia (ET) and myelofibrosis (MF). The majority of patients with CALR mutations have either a 52-base pair deletion (CALR<sup>DEL</sup>) or a 5-base pair insertion (CALR<sup>INS</sup>) in exon 9. This protocol detects CALR<sup>DEL</sup> and CALR<sup>INS</sup> mutations, as well as other less common CALR deletion mutations (e.g. the 61-base pair deletion found in the Marimo cell line). The CALRgenoF1 + CALRgenoR primer set amplifies the muated region of the gene, while the CALRgenoF2 primer is specific to the 5-base pair insertion and, along with CALRgenoR, generates a smaller PCR product. CALR mutations may then be identified via PCR product size discrimination.</p> <p>This protocol is optimized for genomic DNA purified from ACK pellets or cell lines using a commercial kit or our <i>Isolation of Genomic DNA from Peripheral Blood</i> protocol. DNA isolated from colonies grown in methylcellulose may also be used; for best results, pick colonies into 100µl of nuclease-free water and incubate in the thermal cycler at 97°C for 15 minutes with heated lid to lyse cells. Store plates at -20°C.</p> <p>For more information, see Jeong, et al. (2016). Screening PCR versus Sanger sequencing: detection of CALR mutations in patients with thrombocytosis. <i>Annals of Laboratory Medicine</i>, 36(4), 291-299.</p>																				
Materials	<ol style="list-style-type: none"><li>1. Nuclease-free water</li><li>2. Dream<i>Taq</i> green DNA polymerase and 10x PCR buffer</li><li>3. dNTPs (10mM)</li><li>4. Primers (10µM): CALRgenoF1: 5'-GCAGCAGAGAAACAAATGAAGG-3' CALRgenoF2: 5'-GCAGAGGACAATTGTCGGA-3' CALRgenoR: 5'-AGAGTGGAGGAGGGGAACAA-3'</li><li>5. 96-well plates and seals or PCR tube strips</li><li>6. 1.5ml microcentrifuge tube or reagent reservoir</li><li>7. Thermal cycler</li><li>8. 2% DNA gel and low-range DNA ladder (25-700bp)</li></ol>																				
Protocol	A. PCR	Notes																			
1.	<div>Allow DNA samples to thaw on ice. Prepare master mix as follows:<table><tr><th>Reagent</th><th>1 well</th></tr><tr><td>10x PCR buffer</td><td>2 µl</td></tr><tr><td>dNTP mix (10 mM)</td><td>0.4 µl</td></tr><tr><td>CALRgenoF1 primer (10 µM)</td><td>0.5 µl</td></tr><tr><td>CALRgenoF2 primer (10 µM)</td><td>0.5 µl</td></tr><tr><td>CALRgenoR primer (10 µM)</td><td>1 µl</td></tr><tr><td>Nuclease-free H<sub>2</sub>O</td><td>14.44 µl</td></tr><tr><td>DreamTaq DNA polymerase</td><td>0.16 µl</td></tr><tr><td>DNA (~10ng/µl)</td><td>1 µl</td></tr></table></div>	Reagent	1 well	10x PCR buffer	2 µl	dNTP mix (10 mM)	0.4 µl	CALRgenoF1 primer (10 µM)	0.5 µl	CALRgenoF2 primer (10 µM)	0.5 µl	CALRgenoR primer (10 µM)	1 µl	Nuclease-free H <sub>2</sub> O	14.44 µl	DreamTaq DNA polymerase	0.16 µl	DNA (~10ng/µl)	1 µl	<i>Be sure to make enough master mix for at least 1 additional well. For a full 96-well plate, make enough master mix for 110 wells and mix in a reagent reservoir.</i>	
Reagent	1 well																				
10x PCR buffer	2 µl																				
dNTP mix (10 mM)	0.4 µl																				
CALRgenoF1 primer (10 µM)	0.5 µl																				
CALRgenoF2 primer (10 µM)	0.5 µl																				
CALRgenoR primer (10 µM)	1 µl																				
Nuclease-free H <sub>2</sub> O	14.44 µl																				
DreamTaq DNA polymerase	0.16 µl																				
DNA (~10ng/µl)	1 µl																				

2.	Transfer 19µl master mix to each well or tube. Add 1µl DNA per well. Seal plate well and briefly centrifuge in plate spinner.	<i>Always balance plate spinner. For tube strips, spin down using the tube strip adapter for the minicentrifuge.</i>																			
3.	Place plate into thermal cycler and run on the following program: <div style="text-align: center;"> <table border="1"> <tr><td colspan="3">Lid temp: 105°C</td></tr> <tr><td>1x</td><td>94°C</td><td>5 min</td></tr> <tr><td rowspan="3">39x</td><td>94°C</td><td>30 sec</td></tr> <tr><td>64°C</td><td>30 sec</td></tr> <tr><td>72°C</td><td>30 sec</td></tr> <tr><td>1x</td><td>72°C</td><td>7 min</td></tr> <tr><td></td><td>4°C</td><td>∞</td></tr> </table> </div>	Lid temp: 105°C			1x	94°C	5 min	39x	94°C	30 sec	64°C	30 sec	72°C	30 sec	1x	72°C	7 min		4°C	∞	<i>For tube strips, be sure to use the green tube support insert.</i>
Lid temp: 105°C																					
1x	94°C	5 min																			
39x	94°C	30 sec																			
	64°C	30 sec																			
	72°C	30 sec																			
1x	72°C	7 min																			
	4°C	∞																			
4.	Prepare 2% agarose DNA gel(s): <ul style="list-style-type: none"> <li>Prepare 2 large gels for each 96-well plate. Each large gel uses 2 combs and is about 200-250ml. Use 12-25µl EtBr for each gel.</li> <li>For fewer samples, prepare a small gel using the 8-well combs. Each small gel is about 50ml and uses 5µl EtBr.</li> </ul>	<i>Allow the gel to cool somewhat before adding EtBr! Add EtBr directly to the gel tray and use a comb to mix thoroughly.</i>																			
1.	Run PCR and digested products on a 2% agarose gel alongside 8µl of low-range DNA ladder. Load 15-20µl PCR product and/or 25-30µl digested product. <div style="text-align: center;">  </div>	<ul style="list-style-type: none"> <li>CALR<sup>WT</sup> = 357 bp</li> <li>CALR<sup>INS</sup> = 272 bp</li> <li>CALR<sup>DEL</sup> = 302 bp (varies based on patient's specific mutation; this estimate reflects the 52-bp deletion, which is most common)</li> </ul>																			